PHAKIC REFRACTIVE LENS/DEVICE DESCRIPTION:

The Medennium Phakic Refractive Lenses (MPL) are precision made optical devices, available in different models/sizes for intraocular implantation. The optical portion is biconcave or concave-convex and made from silicone elastomer. The haptics are spherical and flexible for ease of insertion, and the overall lens shape is designed to conform to the natural shape of the posterior chamber of the human eye, and to move with the dynamic changes of the eye. Gross centration, if necessary, is achieved via haptic contact with the ciliary sulcus with precise centration achieved through the iris self-centering design and the central optic fenestration (Parachute Technology). The MPL is not fixed to the ciliary sulcus but is allowed by design to float in its self-centered position. Refractive index of the MPL lens is 1.46.

Model	Optic Dia \varnothing_R	Overall Dia \varnothing_T	Overall Width	Diopter D	Increments
100	4.3 to 6.0 mm	10.8 mm	6.0 mm	+1.0 to -28.0	0.5 D
101	4.3 to 6.0 mm	11.3 mm	6.0 mm	+1.0 to -28.0	0.5 D
200	4.5 to 5.8 mm	10.6 mm	6.0 mm	+1.0 to +15.0	0.5 D

INDICATIONS: Visual correction of myopia, hyperopia, and correction of refraction in phakic eyes 18 years of age and older. The lens is intended to be placed in the posterior chamber.

PRECAUTIONS: Do not resterilize the lens by any method. Use only sterile intraocular irrigation solutions and viscoelastics to rinse and/or lubricate the lens. Do not store the lens in direct sun light, or at temperatures over 50°C (122°F). Do not autoclave the lens. Reuse and/or reesterilization may compromise device performance, which could cause serious harm to the patient’s health and safety.

WARNINGS:

1. Physicians considering MPL implantation under any of the following circumstances should weigh the potential risk/benefit ratio. The lens may exacerbate or may interfere with diagnosis or treatment of the conditions below, which may pose an unrelated risk to the patient’s eyesight:
 - a) Recurrent anterior segment or posterior segment inflammation or uveitis.
 - b) Patients in whom the lens may interfere with the ability to observe, diagnose or treat posterior segment disease.
 - c) Surgical difficulties, which might increase the potential for complications (e.g., persistent bleeding, significant iris damage, uncontrolled intraocular pressure or corneal endothelial dystrophy).
 - d) Prior history of retinal detachment, open angle glaucoma or suspect diabetes.
 - e) Patients with only one eye with potentially good vision.
2. One generous peripheral YAG iridotomy (or surgical iridectomy) should be present and confirmed by retroillumination.
3. Inflammation can be an indication of infection and must be treated immediately.
4. As with any surgical procedure, there is risk involved. Potential complications accompanying implant surgery may include, but are not limited to the following: corneal endothelial damage, cystoid macular edema, infection, hypopyon, pupillary block, secondary glaucoma, iris prolapse and uveitis.
5. The long-term effect of MPL implantation has not been determined. Therefore, physicians should continue to monitor MPL patients postoperatively on a regular basis.
6. Physicians considering MPL implantation in patients with macular degeneration should be aware that the patient may not achieve any improvement in central visual acuity following implantation.
7. Patients with ocular pathology (e.g. glaucoma or corneal disease) may not achieve the visual acuity of patients without such problems.
8. Implantation of the MPL in the anterior chamber is unsafe.

CONTRAINDICATIONS:

1. Implantation of the MPL is contraindicated in patients with weak zonula such as previous history of ocular trauma, phacodoneisis, Marfan Syndrome, Weill-Marchesani Syndrome, pseudoexfoliation syndrome, previous history of uveitis and ocular inflammation.
2. Patients with chronic systemic diseases such as Diabetes Melitus, Ulcerative Colitis, Systemic Lupus Erythematosus, Crohn’s Disease, Collagen Vascular disease, Rheumatoid Arthritis, or any bleeding diathesis should be contraindicated.
3. Implantation of MPL when the anterior chamber depth is smaller than 3.0mm when measured from corneal apex to the anterior apex of the crystalline is contraindicated.
4. The lens should not be implanted in eyes that have infection, inflammation, corneal endothelial damage or cystoid macular edema.
5. The safety of MPL implants has not been substantiated in patients with pre-existing ocular conditions (chronic drug miosis, glaucoma, diabetic retinopathy, previous corneal transplant, history of retinal detachment, or iritis, etc.). Physicians considering lens implantation in such patients should explore the use of alternative methods and consider lens implantation only if alternatives are deemed unsatisfactory to meet the need of the patient.

ADVERSE REACTIONS:

As with any surgical procedure, there are risks involved. Based on previous clinical experience, post-operative complications or adverse reactions may occur, including but not limited to the following:

- MPL damage during implantation.
- Traumatic cataract because of crystalline trauma during implantation
- Corneal endothelial damage related to endothelial touch during implantation
- Infection
- Hypopion
- Pupillary block
- Secondary Glaucoma
- Iris prolapse
- Under or over corrections
- MPL decentration and/or subluxation
- Dislocation of the MPL lens into the vitreous cavity
- Glare and /or halos
- Crystalline lens opacification
- Incisional leaking
- Induced astigmatism
- Hyphema
- Hypotony
- Uveitis
- Transient elevated intraocular pressure requiring or not requiring treatment
- Glaucoma
- Ciliary Spasms
- Atonic/mydiatric pupil
- Angle closure
- Cystoid macular edema

DIRECTIONS FOR USE:

1. Administer a retrobulbar block or peribulbar block to reduce eye movement and patient light sensitivity.
2. Examine the label on the lens package for proper lens model, power and expiration date.
3. Open the pouch and remove the lens in a sterile environment.
4. Examine the lens thoroughly. The lens may become electrostatically charged upon opening the package. Examine the lens to ensure that particles are not present on it, rinse if necessary with sterile balanced salt solution.
5. The following folders/injectors are recommended for lens implantation and handling: Dementiev implantation forceps J3100 (e.janach) or 4-2015T/4-2016T (Rumex) or EL44 (Epsilon), Double Ended-Spatula 6-070 (D&K) or Dementiev Spatula J2197A (e.janach) or EP99 (Epsilon), Folding Block EB11 (Epsilon), Dualtec Kit Injector System OD665 (Ophtec), Comport Trimo Ref: T22T7113 RET. Inc.
6. For the purpose of protecting eye structures the following viscoelastics have been clinically validated: Gelbag, Dispersa and Ophthalmic.
7. The surgeon should take care not to traumatize the natural lens or the corneal endothelium to prevent cataract induction and/or endothelial damage.
8. Careful removal of the viscoelastic by aspiration and/or irrigation with balanced salt solution before wound closure to prevent intraocular pressure increase and/or corneal edema related to viscoelastic retention.
9. For the purpose of MPL explantation or repositioning we strongly recommend floating Viscoat on top of the original viscoelastic in order to better protect the endothelium and to use the extraction forceps DK2-796 (D&K)

LENS POWER CALCULATIONS:

The Physician will measure the preoperative manifest & cycloplegic refractions (including cylinder), and the Anterior Chamber Depth (ACD, distance from the vertex of the cornea to the anterior vertex of the natural lens). Calculate the refractive cycloplegic spherical equivalent (CRSE) for the subject taking into account at most 2.0 D of cylinder. Any additional cylinder should be treated by alternative methods since the MPL is designed to compensate for the spherical error only. Report the CRSE and the desired Target Post Operative Refraction (MRSE) to Medennium. Attempt to under-correct, specially if the patient is older. Please note that as the degree of myopia increases the refraction is more difficult to perform. Use of standard contact lenses to knock out most of the correction may be a better method to obtain more accurate refractions.

LENS MODEL SELECTION:

1. MYOPIC: To determine the lens model to use, measure the horizontal “white to white” diameter. If better measurement methods are available (i.e Orbscan) to determine the sulcus to sulcus distance, follow that measurement to fit the MPL so as to prevent artificial vaulting.

Horizontal white to white	MPL Model to implant
≥ 11.3mm	101
< 11.3mm	100

2. HYPEROPIC: All cases should use lens Model 200.

PATIENT REGISTRATION: An ID card is supplied in the lens packaging. This card must be given to the patient with instructions to keep as a permanent record of the implant and to show the card to any eye care practitioner seen in the future.

REPORTING: Adverse reactions and/or potentially sight-threatening complications that may be reasonably regarded as lens related and that were not previously expected in nature, severity or incidence must be reported to Medennium Inc. Irvine, CA USA. Tel: 1-949-789-9000 Fax 1-949-789-9032 E.mail: customerservice@medennium.com. This information is requested from all surgeons in order to document potential long-term effects of intraocular lens implantation.

HOW SUPPLIED/ EXPIRATION DATE: Each lens is supplied sterile, in dry form, in a lens container sealed within a single sterile pouch pack. The lenses are Gamma or Hydrogen Peroxide sterilized. Sterility is guaranteed unless the pouch is damaged or otherwise compromised. Expiration date is indicated on the outside of the lens box. The lens should not be used after the expiration date indicated. Do not store the lens at temperatures below -20°C (-4°F) or above 50°C (122°F).

RETURN/EXCHANGE POLICY: Contact Medennium Inc. for return authorization. Return the lens with proper identification and reason for return. Label return as a biohazard if the lens was in contact with the patient.

DISCLAIMER OF LIABILITY: Medennium Inc. accepts no liability for the choice of method or technique to implant the lens or for the choice of the lens for a particular patient or patient’s condition.

Manufacturer: Medennium Inc. 13 Avenue de Montrouge, E 6, R Normandie, 92340 Bourg La Reine , France

EC REP Dr. Dimitrii Dementiev M.D. 13 Avenue de Montrouge E6. R. Normandie 92340 Bourg la Reine, France.

European Patents: EP0934038, EP1207818, EP1303233, EP1242007
Other worldwide patents pending.

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