RX - Federal (USA) law restricts this device to sale by or on order of a physician.

Instructions for Use: SmartPlug®

Description: Medennium’s SmartPlug® is designed to be inserted into the punctal opening in order to block tear drainage through the canaliculus. Each punctum plug is molded from a thermodynamic acrylic material. The plug is available in one size to fit puncta sizes from 0.4 to 0.8mm.

Anesthesia: Anesthetize the area of the punctum with either a topical anesthetic instilled into the conjunctival sac or an anesthetic soaked, cotton tipped applicator held against the punctum for 30 seconds.

Insertion: SmartPlug® comes preloaded on a disposable silicone carrier. Remove the carrier from the sterile tray and while grasping the top and bottom of the carrier, in a twisting action separate the carrier to expose SmartPlug® see figures 1-4. Grasp the plug with forceps Model 502 and dispose of the carrier. Use the forceps to position the plug across the punctum. Lubrication of the plug with saline solution or artificial tears eases insertion. Gently insert the plug into the punctum 2/3 of its length and remove forceps. Allow a few seconds for plug recovery. This will be obvious by observing the plug move into a more vertical position and retracting into the punctum. Observe and verify proper placement of the plug. The top edge of the plug should be flush with the punctum opening or slightly beneath the level of the punctum opening. Forceps may be used to further push SmartPlug® to obtain the desired position. Always allow enough time for the SmartPlug® to warm up to body temperature before advancing it past 2/3 of its length.

Indications for Use: SmartPlug® is intended for use in patients experiencing dry eye symptoms such as redness, burning, reflex tearing, itching or foreign body sensations which can be relieved by occluding the punctum. SmartPlug® may be used in the treatment of dry eye syndrome and the dry eye components of various ocular surface diseases. When indicated SmartPlug® may be used after ocular surgery to prevent complications due to dry eye and to enhance the retention of topical ocular medications on the eye. Patients experiencing dry eye related contact lens problems may also be aided by SmartPlug®.

Contraindications: SmartPlug® is contraindicated for patients experiencing epiphora, inflammation of the eye lid, and tearing secondary to dacryocystitis with mucopurulent discharge. If the patient experiences irritation, infection, or epiphora after insertion of SmartPlug®, the plug should be removed by flushing the punctum with saline solution or by grasping it with forceps and extracting it from the punctum.

Precautions: Do not re-sterilize SmartPlug® by any method. Use only sterile intraocular irrigation solutions when required. Do not store SmartPlug® in direct sun light, at temperatures over 30°C (86°F) or below -20°C (-4°F).

Care should be used not to perforate the punctum or canaliculus during the insertion of SmartPlug®. Perforation may cause pain, increase risk of infection and cause the plug to fall out more easily due to an enlarged punctal opening or infection. If perforation occurs, delay the placement of the plug until the wound heals.

SmartPlug® may enhance the effect of ocular medication used on the eye. Depending on the type of medications being used, dosage may need to be reduced accordingly. Squeezing of SmartPlug® with sharp-tipped forceps, may cause it to break. Use of rounded forceps or Medennium’s Model 502 is recommended to prevent damage. Be aware that at temperatures above 30°C (86°F) SmartPlug® will recover more rapidly, so it should be implanted as soon as it is removed from its carrier. During implantation avoid touching the device with the lids of the patient to prevent sudden recovery of SmartPlug®.

Reuse and/or re-sterilization may compromise device performance, which could cause serious harm to the patient’s health and safety.

In the presence of infection, treating the infection first will facilitate device removal. The possible risks with use of the plug are irritation, granulomas, infection (Canaliculitis) and/or difficulty in removal. The plug may be removed by irrigation using cannula Model 506 and in exceptional cases by surgical intervention.

How Supplied/Expiration Date: Each SmartPlug® is supplied sterile, in dry form, in a silicone carrier within a single sterile pack. Each outer case contains two Gamma sterilized SmartPlug®. Sterility is guaranteed unless the pack is damaged or otherwise compromised. Expiration date is indicated on the label. SmartPlug® should not be used after the expiration date indicated.

Return/Exchange Policy: Contact Medennium Inc. for return policy. Return the product with proper identification and reason for return. Label return as a biohazard.

Disclaimer of Liability: Medennium Inc. accepts no liability for the choice of method or technique to implant SmartPlug® or for the choice of the product for a particular patient or patient’s condition.

Manufacturer: Medennium Inc., 9 Parker, Suite 150, Irvine, CA 92618, U.S. Tel: (949) 789-9000 Fax: (949) 789-9032. [EC REP], Dr. Dimitrii Dementiev M.D, 13, Avenue de Montrouge, Escalier 6, R Normandie, 92340 Bourg La Reine, France

Clinical Trial Summary: Clinical data were obtained from a prospective, comparative, open-label study of the Medennium, Inc. SmartPlug™ and a silicone control plug (predicate device) for the treatment of dry eye. One hundred and twenty subjects were enrolled. Subjects received a SmartPlug™ into the lower punctum of the randomly selected study eye with the control plug being inserted into the lower punctum of the contralateral eye. Eyes were followed for 90 days after insertion or until the plug was lost or removed (Table 1). One subject was lost to follow-up and one subject exited the study prior to the 90-day visit due to a scheduling conflict.

Table 1. Treatment Failures

Table with 5 columns: Plug lost (fell out) - cumulative SmartPlug Control Plug, 1 Day n/N (%), 7 Days n/N (%), 30 Days n/N (%), 90 Days n/N (%). Rows include SmartPlug Control Plug and Plug removed - cumulative SmartPlug Control Plug.

The primary safety parameter was the rate canaliculitis and dacryocystitis. During the study, these adverse events were not observed in SmartPlug™ or the control plug eyes.

One incidence of each of the following complications was observed in the control group: partially retained silicone plug due to breakage during removal, acute bacterial conjunctivitis, and aggravation of a pre-existing inflammatory lesion of the lower lid. One incidence of each of the following complications was observed in SmartPlug™ group: chemosis and throbbing pain. Epiphora was reported in both eyes of two subjects. Allergic conjunctivitis was reported in both eyes of one subject. Transitory monocular diplopia was reported in both eyes of one subject.

Effectiveness was evaluated by Rose Bengal staining (Table 2). The surface of the eye was divided into three zones, and staining in each zone was graded on a scale from 0 to 3 (none to severe). For each eye, the sum of these three scores at each follow-up visit was compared to baseline. The difference in improvement between the two eyes of each subject was taken into account in the statistical equivalence test for non-inferiority.

Additional effectiveness parameters were similarly analyzed (Tables 3 and 4). Basal tear secretion was evaluated with anesthesia and measured as millimeters of wetting of Schirmer strips after five minutes. Fluorescein staining scores were the sum of grades on a scale from 0 to 3 (none to severe) in five ocular surface zones.

Table 2. Rose Bengal Staining: Paired Eye Analysis

Table with columns: Device, Baseline (N, Mean S.D.), 1 Day (N, Mean S.D.), 7 Days (N, Mean S.D.), 30 Days (N, Mean S.D.), 90 Days (N, Mean S.D.). Rows include Smart Plug, Control Plug, Equivalence Test P-Value, and Equivalent.

Table 3. Basal Tear Secretion: Paired Eye Analysis

Table with columns: Device, Baseline (N, Mean S.D.), 1 Day (N, Mean S.D.), 7 Days (N, Mean S.D.), 30 Days (N, Mean S.D.), 90 Days (N, Mean S.D.). Rows include Smart Plug, Control Plug, Equivalence Test P-Value, and Equivalent.

Table 4. Fluorescein Staining: Paired Eye Analysis

Table with columns: Device, Baseline (N, Mean S.D.), 1 Day (N, Mean S.D.), 7 Days (N, Mean S.D.), 30 Days (N, Mean S.D.), 90 Days (N, Mean S.D.). Rows include Smart Plug, Control Plug, Equivalence Test P-Value, and Equivalent.

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