

SmartPLUG® A simple solution for dry eye

SmartPLUG® gently conforms to the contours of the lacrimal canal with no protruding parts to cause irritation or discomfort

made from an exclusive thermodynamic, hydrophobic acrylic polymer • transforms from a rigid rod into a soft gel-like plug at body temperature • aesthetically pleasing and more comfortable

SmartPLUG® is a thin 6 mm acrylic rod that fits easily into the punctum at room temperature. Topical anesthesia is optional.

Once in contact with the body, SmartPLUG® shrinks in length and expands in width, reshaping itself into a soft, comfortable, gel-like plug.

SmartPLUG® accommodates to the size and shape of most lacrimal canals, eliminating the need for measurement. Dilation is optional.

One in place, nothing protrudes from the punctum, making SmartPLUG® more aesthetically pleasing and comfortable.

If desired, SmartPLUG® can be removed by irrigation.

SmartPLUG® before



SmartPLUG® after



www.medennium.com
Toll Free: 1-888-727-6100


Medennium
Think Better & See Different

Medennium
9 Parker Suite 150
Irvine, CA 92618 USA