


Smart Light Transition

One-Piece Foldable Aspheric IOL

MATRIX
Acrylic® **AURIUM**



This product is CE marked and is only available in the European Union and jurisdictions that recognize CE marking.

Parameters

Material	Hydrophobic acrylic
Power	0.0D to +34.0D
Haptic	L-shape
A-Constant/IOLMaster	118.4/119.0
Optical Diameter	6mm
Overall Diameter	13mm

Smart Light Transition	
Clear-to-yellow	10 seconds
Yellow-to-clear	30 seconds
Incision size	2.4-2.8mm
Refractive index	1.56(20°C)
Spherical Aberration	-0.1 μm



This product is CE marked and is only available in the European Union and jurisdictions that recognize CE marking.



One-Piece Foldable Aspheric IOL

403

Medennium

CE


This product is CE marked and is only available in the European Union and jurisdictions that recognize CE marking.

One-Piece Foldable Aspheric IOL Model 403



Matrix Acrylic Lens may be implanted
with an injector through a 2.4-2.8mm incision

Parameters



Optic: biconvex aspheric

A-Constant: 118.4

Material: Hydrophobic acrylic

IOLMaster: 119.0

Haptic: L-shape

Refractive index: 1.56(20°C)

Power: 0.0D to +34.0D

Spherical Aberration: -0.1 μm


This product is CE marked and is only available in the European Union and jurisdictions that recognize CE marking.

Medennium

www.medennium.com



This product is CE marked and is only available in the European Union and jurisdictions that recognize CE marking.



Three-Piece Foldable Spheric IOL

401


Medennium

CE


This product is CE marked and is only available in the European Union and jurisdictions that recognize CE marking.

Three-Piece Foldable Spheric IOL Model 401



Matrix Acrylic Lens may be implanted
with an injector through a 2.8-3.2mm incision

Parameters



Optic: Fully symmetric biconvex

A-Constant: 118.5

Material: Hydrophobic acrylic

IOLMaster: 119.2

Haptic: PVDF-Modified C

Refractive index: 1.56(20°C)

Power: 0.0D to +30.0D

Medennium

www.medennium.com





This product is CE marked and is only available in the European Union and jurisdictions that recognize CE marking.

Smart Light Transition

Three-Piece Foldable Spheric IOL

MATRIX
Acrylic® **AURIUM**




This product is CE marked and is only available in the European Union and jurisdictions that recognize CE marking.



Parameters

Material	Hydrophobic acrylic
Power	0.0D to +30.0D
Haptic	PVDF-Modified C
A-Constant/IOLMaster	118.5/119.2
Optical Diameter	6mm
Overall Diameter	12.5mm

Smart Light Transition	
Clear-to-yellow	10 seconds
Yellow-to-clear	30 seconds
Incision size	2.8-3.2mm
Refractive index	1.56(20°C)



This product is CE marked and is only available in the European Union and jurisdictions that recognize CE marking.